

## Application note

# Oligomerization of Interleukin-2

A commercially available human interleukin-2 is assessed by microfluidic diffusional sizing across a dilution series. The hydrodynamic radius is observed to increase with increasing concentration, in a way which suggests a monomer-trimer equilibrium with positive cooperativity is established.

### Authors

Maren Butz<sup>1</sup>, Sean Devenish<sup>1</sup>

[www.fluidic.com](http://www.fluidic.com)



## Introduction

Interleukin-2 (IL-2) is a lymphokine with a range of crucial roles in immune activation and homeostasis (1). It is produced by a subset of activated T-lymphocytes and acts as a growth factor for T cells and NK cells (2), and promotes proliferation, activation-induced cell death, antibody secretion, and cytokine production in various cell types (3). Intravenous treatment with high dose IL-2 has been used in patients with certain cancers (4) (5).

IL-2 has a molecular weight of 15.5 kDa and is globular, which corresponds to a size of 2.04 nm for a monomeric form (6).

Here we use Microfluidic Diffusional Sizing (MDS) on a Fluidity One instrument to assess the hydrodynamic radius ( $R_h$ ) over a dilution series. The tests use a commercially available lyophilized form with no additives. We observe that the  $R_h$  increases with increasing concentration, in a manner which indicates a monomer-trimer equilibrium with positive cooperativity.

This result challenges much of the modern literature which assumes a monomeric form (7) (8) (9), or in some instances a dimer (10). One study by Siddarth et al uses analytical ultracentrifugation and produced data which cannot rule out a monomer-trimer equilibrium, though a monomer-dimer model is eventually concluded (11). It should be noted that a different IL-2 source and buffer is used in this study.

## Methods

50  $\mu\text{g}$  of lyophilized human IL-2 recombinantly expressed in *E. coli* (Sigma-Aldrich code SRP3085) was dissolved in 10 mM acetic acid at 1 mg/mL.

A 10  $\mu\text{L}$  aliquot of this stock was diluted with 23.3  $\mu\text{L}$  of 138.6 mM acetic acid to produce a solution of nominally 300  $\mu\text{g}/\text{mL}$  in 100 mM acetic acid. (The subsequent measured concentrations of the dilution series being consistently higher than expected suggests that this initial stock was > 300  $\mu\text{g}/\text{mL}$ , likely due to the purchased 50  $\mu\text{g}$  containing more than expected.)

This 300  $\mu\text{g}/\text{mL}$  solution was diluted in 2-fold dilution series using 100 mM acetic acid to give 7 varying concentrations, with a lower limit of 4.7  $\mu\text{g}/\text{mL}$ .

Each concentration was analysed using a Fluidity One prototype instrument to measure the  $R_h$ . Tests were performed in triplicate, with the values reported below being an average of the repeats.

The observed  $R_h$  change over changing concentration was then assessed against different models to infer the nature of the oligomerisation.

## Results

The concentrations of the dilution series along with the measured  $R_h$  and concentration are shown in Table 1.

Given the molecular weight of IL-2, the theoretical  $R_h$  for a monomer is 2.04 nm, dimer is 2.62 nm and trimer is 3.03 nm. As such the sizes observed at higher concentrations are consistent with trimer formation.

Figure 1 shows fitting of the data points to a monomer-dimer equilibrium model (see appendix 1), which does not accurately represent the data.

Table 1: Results of MDS analysis on the IL-2 dilution series

Nominal Concentration (µg/mL)	Average Measured $R_h$ (nm)	Average Measured Concentration (µg/mL)	Notes
4.7	1.87	4.32	-
9.4	2.00	10.75	Average of 2 runs
18.8	2.29	22.53	-
37.5	2.85	40.47	-
75	2.98	88.90	-
150	2.91	190	Average of 2 runs
300	3.05	411	-




Figure 1: Results of dilution series analysis fitted to both a monomer-dimer equilibrium model, and a Hill Equation

It is clear that the monomer-dimer fitting fails at capturing the steepness of the curve, indicating that there may be cooperative behaviour.

Fitting the data instead to a Hill equation provides a more representative fit, as shown in Figure 1. Akaike's Information Criterion Test also shows that the Hill equation is ~3000-fold more likely than the monomer-dimer equilibrium model.

The steepness of the slope observed in this fitting indicates there is positive cooperativity in the oligomerisation, as indicated in Table 2.

Table 2: Hill equation fitting of the dilution series data  
 No cooperativity;  $n = 1$   
 Positive cooperativity;  $n > 1$

$y = MR_h + (DR_h - MR_h) * x^n / (K_d^n + x^n)$	Result
$MR_h$ (nm)	$1.88 \pm 0.05$
$DR_h$ (nm)	$3.00 \pm 0.04$
$K_d$ (µg/mL)	$22 \pm 2$
$K_d$ (µM)	$1.42 \pm 0.13$
$n$	$3.3 \pm 0.7$
$R^2$	0.961

## Conclusion

A change in the  $R_h$  of IL-2 on dilution has been observed using microfluidic diffusional sizing on a Fluidity One instrument. The  $R_h$  changed from 3.1 nm at 300 µg/mL concentration, to 1.9 nm at 4.7 µg/mL.

The change in size indicates that oligomerization occurs at concentrations above 20 µg/mL. The size observed at 300 µg/mL is consistent with a trimer.

When the data is analysed, a Hill equation provides a better fit compared to a monomer-dimer equilibrium model. The steepness of the curve on the Hill equation fit indicates positive cooperativity in the assembly of the trimers.

## References

1. Gaffen, S. L. and Liu, K. D., 2004. Overview of interleukin-2 function, production and clinical applications, *Cytokine*. 28(3). p.109.
2. Schorle, H., Holtschke, T., Hünig, T., Schimpl, A. and Horak, I., 1991. Development and function of T cells in mice rendered interleukin-2 deficient by gene targeting, *Nature*. 15. p.621.
3. Liao, W., Lin, J. and Leonard, W., 2012. IL-2 Family Cytokines: New Insights into the Complex Roles of IL-2 as a Broad Regulator of T helper Cell Differentiation, *Current Opinions in Immunology*. 23(5). p.598.
4. Atkins, M. B., Lotze, M. T., Dutcher, J. P., Fisher, R. I., Weiss, G., Margolin, K., Abrams, J., Sznol, M., Parkinson, D., Hawkins, M., Paradise, C., Kunkel, L., Rosenberg, S. A., 1999. High-dose recombinant interleukin 2 therapy for patients with metastatic melanoma: analysis of 270 patients treated between 1985 and 1993, *Journal of Clinical Oncology*. 17(7). p.2105.
5. Fyfe, G., Fisher, R. I., Rosenberg, S. A., Sznol, M., Parkinson, D. R. and Louie, A. C., 1995. Results of treatment of 255 patients with metastatic renal cell carcinoma who received high-dose recombinant interleukin-2 therapy, *Journal of Clinical Oncology*. 13(3). p.688.
6. Fluidic Analytics Ltd. How do you convert hydrodynamic (Stokes) radius to molecular weight? [Online] April 2018. <https://www.fluidic.com/resources/faq/convert-hydrodynamic-radius-to-mw/>
7. Wrenshall, L. E., Clabaugh, S. E., Cool, D. R., Arumugam, P., Grunwald, W. C., Smith, D. R., Liu, G. C. and Miller, J. D., 2014. Identification of a Cytotoxic Form of Dimeric Interleukin-2 in Murine Tissues, *PLoS One*. 9(7).
8. Liang, S. M., Allet, B., Rose, K., Hirschi, M., Liang, C. M. and Thatcher, D. R., 1985. Characterization of human interleukin 2 derived from *Escherichia coli*, *Journal of Biochemistry*. 229(2). p.429.
9. Ricci, M. S., Sarkar, C. A., Fallon, E. M., Lauffenburger, D. A. and Brems, D. N., 2003. pH Dependence of structural stability of interleukin-2 and granulocyte colony-stimulating factor, *Protein Science*. 12(5). p.1030.
10. Kaplan, D. R., 1994. Delivery of interleukin 2 for immunotherapy, *Journal of chromatography B: Biomedical Sciences and Applications*. 662(2). p.315.
11. Advant, S. J., Braswell, E. H., Kumar, C. V. and Kalonia, D. S., 1995. The effect of pH and temperature on the self-association of recombinant human interleukin-2 as studied by equilibrium sedimentation, *Pharmaceutical Research*. 12(5). p.637.

## Appendix 1 - Data fitting equations

Fitting equation for Monomer-Dimer equilibrium:

$$K_d = [M]^2 / [D]$$

Total concentration of IL-2:

$$x$$

Concentration of IL-2 monomer:

$$[M] = x - 2z$$

Concentration of IL-2 dimer:

$$[D] = z$$

Quadratic equation:

$$z = \frac{(4x + K_d) - \sqrt{(4x + K_d)^2 - 16 \times x^2}}{8}$$

$$\text{Equation for Average } R_h: \text{ Ave}R_h = \left( \frac{(4x + K_d) - \sqrt{(4x + K_d)^2 - 16 \times x^2}}{8} \right) \frac{2}{x} \times DR_h + \left( 1 - \left( \frac{(4x + K_d) - \sqrt{(4x + K_d)^2 - 16 \times x^2}}{8} \right) \frac{2}{x} \right) \times MR_h$$

Average  $R_h$ :

$\text{Ave}R_h$

Monomer  $R_h$ :

$MR_h$

Dimer  $R_h$ :

$DR_h$